

Santa Rosa Junior College
Radiologic Technology Program

C-Arm Orientation Checklist

Student Name: _____	Semester: _____
---------------------	-----------------

Locate and/or operate	completed	N/A
- Brakes and steering mechanisms.		
- Connect C-Arm unit to monitors.		
- ON/OFF switch/button.		
- Exposure technique control buttons/knobs.		
- Low dose and boost control.		
- Contrast and brightness control on monitor.		
- Collimation control.		
- Image orientation control.		
- Fluoro timer reset.		
- Movement control levers/handles.		
- Image save/store buttons.		
- Exposure switches (hand, foot), controls.		
- Hard copy devices.		
- Data entry using keyboard.		
- Annotate data before and after procedure.		
- Storage location.		
- Send images to PACS		
Radiation Protection		
- Understands how surgical cases are ordered.		
- Only expose when ordered by the physician.		
- Make sure all personnel are wearing protective aprons.		
Advanced Procedures (if applicable)		
- Cine radiography		
- Road mapping		
- Image subtraction		
- Peak opacification		
- Storing of images and cine		

Comments:

Evaluating R.T. _____	Date: _____
-----------------------	-------------